


# ANTI BEDSORES MATTRESS NINFA



## MAIN SPECIFICATIONS

- Mattress developed for the prevention and treatment of decubitus sores;
- Mattress articulated and not reversible;
- Its constitution allows to minimize the pressure on the skin to avoid scarring and accelerate the healing of existing;
- The mattress Ninfa ensures ergonomically correct body posture during sleep;
- There are two versions mattress Ninfa:
  - Ninfa 12: 7cm soyafoam with the support base of over 5 cm Memory Foam in the surface layer;
  - Ninfa 14: 9cm soyafoam the support base of over 5 cm Memory Foam in the surface layer;
- Outer shell: waterproof bielastic to liquids and gas-permeable, breathable, allowing the breath of the interior, facilitating the hygiene of the mattress. Fire-resistant according to BS 7175 1989, Crib 5. The tissue can be washed to a temperature of 95 ° C with neutral detergents;
- Memory Foam: viscoelastic foam, extremely moldable and cozy, allowing a proven reducing pressure on the skin and blood vessels, contributing to better blood circulation. The viscoelastic foam adapts fully to the body shape, providing a sense of lightness that translates into a very restful sleep;
- Soyafoam: foam used as support base. It is an innovative Polyurethane, because it is an environmentally friendly foam made based on soybean;

- Dimensions:
  - Ninfa 12: 190x90x12cm;
  - Ninfa 14: 190x90x14cm.
- Notes:
  - References mention mattresses for single beds (up to 190 x 90 cm);
  - There is also for couples (up to 200 x 150 cm);
  - They make up mattresses for measures (by budget).

Code	Description	Selling unit
30009014	ANTI BEDSORES MATTRESS NINFA	CX. 1