

## ANTI BEDSORES SEAT BLUE GEL



### MAIN SPECIFICATIONS

- The BLUE GEL seat helps to distribute the pressure on the skin, decreasing the risk of developing pressure ulcers, since they move with the users, not against them;
- Made of polymer gel, does not drip, does not contain latex, plasticizers or silicones, common causes of allergic reactions;
- Washable and disinfected;
- Suitable for hospital environment;
- They have washable covers and disinfected;
- The top in velor, pleasant to the touch and anti-slip base;
- Size available:
  - 40x41x1,8cm;
  - 43x41x1,8cm;
  - 46x41x1,8cm;
  - 46x41x2,8cm;
  - 50x41x2,8cm.

Code	Description	Selling unit
------	-------------	--------------

