

## MATTRESS HERA MEMORY



### MAIN SPECIFICATIONS

- Anti bedsores mattress Hera Memory with 14cm tall;
- Base in polyurethane foam with a density of 28kg / m<sup>3</sup>, high 9cm and top with viscoelastic foam 5cm, with density 50kg / m<sup>3</sup>, totaling 14cm tall;
- Bielastic waterproof cover to liquids and gas-permeable, breathable, allows breathing of the interior facilitates the hygiene of the mattress, fireproof protection;
- Also available Hera Memory with 12cm tall (+ 8cm 4cm viscoelastic polyurethane);
- Suitable for use on articulated beds to prevent bedsores;
- 2 year warranty against manufacturing defects, except for rapid wear of components;
- Estimated Delivery time: 6-8 weeks;
- Two models are available:
  - Hera Mem 14: 190x90x14cm;
  - Hera Mem 12: 190x90x12cm.

Code	Description	Selling unit
30009021	MATTRESS HERA MEMORY 190X90X12 CM	CX. 1
30009022	MATTRESS HERA MEMORY 190X90X14 CM	CX. 1

FACTORMED  
Produtos Médicos & Hospitalares