

# ALTERNATE PRESSURES PAD COMBO 100 PLUS



## MAIN SPECIFICATIONS

- Combo plus 100 is an alternating pressure air cell cushion for the prevention of pressure ulcers;
- Ideal for people with reduced mobility and for those who use a wheelchair to move;
- System composed with four air cells for a foam base;
- Structure which greatly relieves pressure when sitting on a hard surface;
- This pad can be used connected to the electrical current, as with batteries, because the NI-MH batteries are rechargeable and can a battery life of up to 8 hours recharged to the maximum;
- Allows the user to use with complete peace of mind while you are away from home;
- Preventive treatment of pressure ulcers;
- It supports a maximum weight of 120kg;
- Dimensions: 43x43x10cm;
- Weight: 960g.

Code	Description	Selling unit
63 AS AE EMA	ALTERNATE PRESSURES PAD COMBO 100 PLUS	CX. 1

FACTORMED  
Produtos Médicos & Hospitalares