


SOFTSEAT MEMORY FOAM



MAIN SPECIFICATIONS

- The SOFTSEAT in Memory Foam anatomical allow patients to spend long periods sitting with reduced risk of bedsores, assisting in much their comfort;
- May be used as wheelchairs seat, but also for ordinary chairs, for users with problems in positioning or poor flexibility in the joints;
- Ideal for improving posture and positioning;
- Indicated for pain relief in sacrocóccigea and perineal area (postpartum, postoperative pelvic surgeries and anorectal) hemoroidas, bed sore prevention in the sacrum coccyx area;
- Lightweight and easy to carry, available in various thicknesses and materials;
- Cover with three finishes to choice: polisoft (black), tough and looking skin for comfort; bielastic blue, liquid impermeable and breathable, washable and disinfected; 3D bielastic with top with an additional fabric layer to 3D, which avoids the undesirable effects of perspiration;
- Interior memory foam (MF) molded, anatomical;
- Bedsores prevention models:
 - SOFTSEAT comfort black cape - polisoft: 42x42x7cm;
 - SOFTSEAT anti-bed sore blue cover - bielastic: 42x42x7cm;
 - SOFTSEAT anti-bed sore blue cover - bielastic w/ 3D: 42x42x7cm; 52x42x10cm;

- Ideal models for localized bedsores:
 - SOFTSEAT anti-bedsores coccyx- bielastic cover w/ 3D: 42x42x10cm;
 - SOFTSEAT anti-bedsores horseshoe - bielastic cover w/ 3D: 42x42x10cm;
 - Center hole anti-bedsores SOFTSEAT - bielastic cover w/ 3D: 42x42x10cm;
- Kits for chairs:
 - Blue 3D anatomical lumbar support (also available in polisoft);
 - Anti-bedsores kit for wheelchairs bielastic cover.

Code	Description	Selling unit
63 AS CF PO	SOFTSEAT COMFORT BLACK CAPE- POLISOFT 42X42X7 CM	CX. 1
63 AS AE BI 07 PU	SOFTSEAT ANTI-BEDSORES BLUE CAPE- BIELASTIC 42X42X7 CM	CX. 1
63 AS AE BI 07	SOFTSEAT ANTI- BEDSORES BLUE CAPE- BIELASTIC W/ 3D 42X42X7CM	CX. 1
63 AS AE BI 10	SOFTSEAT ANTI- BEDSORES BLUE CAPE- BIELASTIC W/ 3D 52X42X10 CM	CX. 1
63 AS AE CX	SOFTSEAT ANTI- BEDSORES COCCYX- BIELASTIC CAPE W/ 3D 42X42X10 CM	CX. 1
63 AS AE FR	SOFTSEAT ANTI- BEDSORES HORSESHOE SHAPE- BIELASTIC CAPE W/ 42X42X10 CM	CX. 1
63 AS AE BR	SOFTSEAT ANTI- BEDSORES W/ CENTRAL HOLE- BIELASTIC CAPE W/ 3D 42X42X10 CM	CX. 1
63 AP LB 3D	LUMBAR ANATOMICAL SUPPORT 3D BLUE	CX. 1
63 AP LB PO	LUMBAR ANATOMICAL SUPPORT POLISOFT	CX. 1
63 KT AE BI	ANTI- BEDSORES KIT FOR WHEELCHAIRS- BIELASTIC CAPE (ALSO IN 3D)	CX. 1