

# HORSESHOE SHAPE SEAT BLUE GEL



## MAIN SPECIFICATIONS

- The BLUE GEL seat helps to distribute the pressure on the skin, decreasing the risk of pressure sores because it moves with the users, not against them;
- Made of polymer gel, does not drip, does not contain latex, plasticizers or silicones, common causes of allergic reactions;
- Washable and disinfected;
- Suitable for hospital environment;
- The seats also have washable covers and disinfected;
- The top in velor, pleasant to the touch and anti-slip base.

Code	Description	Selling unit
63 AS 4FR 1540	HORSESHOE SHAPE SEAT BLUE GEL 40X41X1,5CM	CX. 1
63 AS 4FR 1543	HORSESHOE SHAPE SEAT BLUE GEL 43X41X1,5 CM	CX. 1
63 AS 4FR 1546	HORSESHOE SHAPE SEAT BLUE GEL 46X41X1,5 CM	CX. 1
63 AS 4FR 1550	HORSESHOE SHAPE SEAT BLUE GEL 50X43X1,5 CM	CX. 1