

COMPLETE BODY PAD BLUE GEL



MAIN SPECIFICATIONS

- The BLUE GEL seat helps to distribute the pressure on the skin, decreasing the risk of pressure sores because it moves with the users, not against them;
- Made of polymer gel, does not drip, does not contain latex, plasticizers or silicones, common causes of allergic reactions;
- Washable and disinfected;
- Suitable for hospital environment;
- The seats also have washable covers and disinfected;
- The top in velor, pleasant to the touch and anti-slip base.
- Dimensions:
 - Head: 50x23x1,5cm;
 - Body: 100x50x1,5cm;
 - Lower limbs: 50x50x1,6cm.

Code	Description	Selling unit
63 AP AE GE CC	COMPLETE BODY PAD BLUE GEL	CX. 1