


## BLUE GEL SEATS



### MAIN SPECIFICATIONS

- The BLUE GEL seat helps to distribute the pressure on the skin, decreasing the risk of pressure sores because it moves with the users, not against them;
- Made of polymer gel, does not drip, does not contain latex, plasticizers or silicones, common causes of allergic reactions;
- Washable and disinfected;
- Suitable for hospital environment;
- The seats also have washable covers and disinfected;
- The top in velor, pleasant to the touch and anti-slip base.
- Models:
  - BLUE GEL for the face: 25x20x10cm;
  - BLUE GEL: occipital (20x7x5cm); pediatric (14x5x3cm; 14x7x2,5cm);
  - BLUE GEL horseshoe shape: occipital (20x7x5cm); pediatric (10x4x2,5cm);
  - BLUE GEL universal head: 25x23x6cm; 25x18x3cm;
  - BLUE GEL oftalmic head: 23x23x4cm; 24x22x7cm;
  - BLUE GEL cervical head: 18x18x8cm;
  -

- BLUE GEL body: 50x20x7cm;
- BLUE GEL wedge: 25x15x6cm;
  - BLUE GEL semi- cylindrical: 30x15x7cm; 30x10x4cm; 50x15x7cm;
  - BLUE GEL arm/ leg: 40x15x4cm;
  - BLUE GEL heel: 18x10x7cm.

Code	Description	Selling unit
63 AC AE GE FC2	BLUE GEL FOR THE FACE 25X20X10 CM	CX. 1
63 AC AE GE OT 2	BLUE GEL OCCIPITAL 20X7X5 CM	CX. 1
63 AC AE GE OT 3	BLUE GEL PEDIATRIC 14X5X3 CM	CX. 1
63 AC AE GE OT 7	BLUE GEL PEDIATRIC 14X7X2,5 CM	CX. 1
63 AC AE GE OT F2	BLUE GEL HORSESHOE SHAPE OCCIPITAL 20X7X5 CM	CX. 1
63 AC AE GE OT F4	BLUE GEL HORSESHOE SHAPE PEDIATRIC 10X4X2,5 CM	CX. 1
63 AC AE GE UN 1	BLUE GEL UNIVERSAL HEAD 25X23X6 CM	CX. 1
63 AC AE GE UN 2	BLUE GEL UNIVERSAL HEAD 25X18X3 CM	CX. 1
63 AC AE GE OF 1	BLUE GEL OFTALMIC HEAD 23X23X4 CM	CX. 1
63 AC AE GE OF 3	BLUE GEL OFTALMIC HEAD 24X22X7 CM	CX. 1
63 AC AE GE CV	BLUE GEL CERVICAL HEAD 18X18X8 CM	CX. 1
63 AT AE GE	BLUE GEL BODY 50X20X7 CM	CX. 1
63 AC AE GE 2	BLUE GEL WEDGE 25X15X6 CM	CX. 1
63 AI AE GE 1	BLUE GEL SEMI- CYLINDRICAL 30X15X7 CM	CX. 1
63 AI AE GE 2	BLUE GEL SEMI- CYLINDRICAL 30X10X4 CM	CX. 1
63 AI AE GE 8	BLUE GEL SEMI- CYLINDRICAL 50X15X7 CM	CX. 1
63 AP AE GE PB 2	BLUE GEL ARM/ LEG 40X15X4 CM	CX. 1
63 AP AE GE CL 1	BLUE GEL HEEL 18X10X7 CM	CX. 1