

LIBERTY SEAT



MAIN SPECIFICATIONS

- Waterproof or velvet cover for comfort, washable;
- Interior with air cells surrounded by gel and air. Adjustable each user through an air valve;
- Splitting the seat into two compartments decreases greatly the continuous tendency to slide forward;
- The front compartment acts as a wedge to inhibit the occurrence of shear forces.
- Accessories: waterproof cape (for incontinence) or standard cover (confair).

Code	Description	Selling unit
63 AS AE PT 40/45	LIBERTY SEAT	CX. 1