

PELVIC STABILIZER LIBERTY



MAIN SPECIFICATIONS

- The pelvic stabilizer is a hard foam base that is placed in the bones (ischial tuberosities) on the seat;
- This raises the seat surface under the thighs, thereby reducing the tendency of the pelvis to rotate backwards (posteriorly positioned), lowering the shear forces and pressure on the sacrum;
- Shall be used on the seat LIBERTY in combination, provide a suitable support.

Code	Description	Selling unit
99 63 E ET PV 40/45	PELVIC STABILIZER LIBERTY	CX. 1