

# ACADEMY ADJUSTER SEAT



## MAIN SPECIFICATIONS

- Seat for symmetric and asymmetric people, even exceptional cases;
- With waterproof or comfort (optional) cape and inside consisting of several compartments with air cells to ensure the most appropriate placement of the user;
- It offers the possibility to create an individual solution for users with functional body asymmetries and/ or severe anatomical;
- The front pad consists in two compartments, which can be adjusted in density;
- Recommended as a good solution to the hemiplegic patients.

Code	Description	Selling unit
63 AS AD 40/45	ACADEMY ADJUSTER SEAT	CX. 1