


MATTRESS URANO



MAIN SPECIFICATIONS

- Mattress in polyurethane with a density of 28kg / m³;
- The setting in jagged or wavy cubes gives is a preventive valence bedsores by decreasing the pressure and for the air flow through its channels and by simple body movements, adapting anatomically;
- Outer layer:
 - The mattress comes with a waterproof bielastic cover to liquids and gas-permeable, breathable that allows the breath of the interior, facilitating the hygiene of the mattress;
 - Fire-resistant in accordance with BS 7175 1989, Crib 5.
 - The tissue can be washed to a temperature of 95 ° C with neutral detergents;
- Dimensions: 190x90x14cm.

Code	Description	Selling unit
61 CL HP UR CB	MATTRESS URANO	CX. 1