

MATTRESS DOMUS IV



MAIN SPECIFICATIONS

- The DOMUS IV is a bedsore prevention system with alternating pressure digitally controlled for hospital and home use;
- Sensory control technology IV DOMUS pressure ensures a consistent and continuous support of low pressure and to prevent unwanted incidents emptying nor recommended for the prevention of bedsores;
- The DOMUS IV offers with a high quality, cost effectiveness, as well as a variety of functions for easy and simple use;
- Compressor:
 - Voltage: AC 100-120V / 50Hz; AC 220-240V / 60Hz;
 - Air flow: 8 liters per minute;
 - Pressure range: 15 - 45 mmHg;
 - Cycle time: 10, 15, 20, 25 selectable;
 - Dimensions (W x L x H): 29 x 20 x 11cm;
 - Weight: 3,7kg;
- Mattress:
 - Structure: 20 cells 30cm alt.;
 - Material: cover (nylon / PU dark blue); cells (nylon / PU light blue); base (nylon / dark gray PVC);

- Dimensions (W x L x H): 200 x 90 x 25.4cm;
- Maximum weight per patient: 200kg.

| Code | Description | Selling unit |
|----------------|-------------------|--------------|
| 61 CL AE DM IV | MATTRESS DOMUS IV | CX. 1 |