

MATTRESS DOUBLE LATEX



MAIN SPECIFICATIONS

- Hospital mattress constructed 100% of latex, firm and resilient material which allows a correct support of the column and facilitates blood circulation;
- The latex is a natural product with very porous optimal ventilation latex being a hypoallergenic product that naturally prevents the growth of bacteria, fungi and mites therein;
- The mattress DOUBLE LATEX still has a surface layer of 2cm Climalátex, which is a very soft latex further ventilated radicional the latex, facilitating the blood circulation;
- It has a great capacadadede ventilation, allowing the quick evaporation of the perspiration, remaining dry and moisture.
- Outer layer:
 - Case in Duvalcare tissue with a composition of 50% Cotton, 25% Polyester 25% Nylon filler with Duvalcare treatment;
 - The Duvalcare fabric has anti-bacterial, anti-fungal and anti-mite;
 - Its characteristics do not change after washing because the wire treatment is permanent;
 - Smooth and supple texture with natural ability to sweat and easy to wash;
 - The fabric is breathable Duvalcare which prevents the occurrence of moisture and odors;
 - It can be washed to a maximum temperature of 60 °C in delicate wash cycle; It can also be dry-

cleaned;

- You can choose between three densities depending on body weight or user preferences:
 - Normal: Latex with a density of 60kg/ m³ for users weighing up to 70kg or preference for low consistency mattresses. Climalatex with density 65,0kg/ m;
 - Firm: Latex with a density of 70kg/ m³, for users with a body weight up to 85kg, or preference for medium consistency mattresses. Climalatex with density 65,0kg/ m.
 - Extra firm: Latex with a density of 80kg/ m³, for users with body weight above 85kg, or preference for high consistency mattresses. Climalatex with density 65,0kg/ m.
- Dimensions: 190x90x15/ 18cm;
- Notes:
 - References mention mattresses for single beds (up to 190 x 90 cm);
 - There is also for couples (up to 200 x 150 cm);
 - We accept orders of mattresses with specific measures (by budget).

Code	Description	Selling unit
61 CL CF DL15 N	MATTRESS DOUBLE LÁTEX 190X90X15 CM	CX. 1
61 CL CF DL18 N	MATTRESS DOUBLE LÁTEX 190X90X18 CM	CX. 1