

# MATTRESS SOFT LIGHT AND SOFT VENUS



## MAIN SPECIFICATIONS

- Soft light: comfort mattress built with a surface layer of 2 cm and a core Memory Foam of 12 or 15cm latex;
- Soft Venus: comfort mattress built with a surface layer of Memory Foam 5cm and a core 12 of latex;
- The surface layer of Memory Foam is a viscoelastic foam, extremely mouldable and cozy, allowing a proven reducing pressure on the skin and blood vessels, contributing to better blood circulation. The viscoelasticity adapts fully to the body shape, providing a sense of lightness that translates into a very restful sleep. Latex core is a natural material, porous, which provides an excellent ventilation;
- It is a hypoallergenic material that naturally prevents the growth of bacteria, fungi and mites therein;
- It is an elastic material, allowing a correct support of the column;
- It is elastic and flexible, facilitating the blood circulation;
- Outer layer:
  - Case in Duvalcare tissue with a composition of 50% Cotton, 25% Polyester 25% Nylon filler with Duvalcare treatment;
  - The Duvalcare fabric has anti-bacterial, anti-fungal and anti-mite;
  - These characteristics do not change after washing because the wire treatment is permanent;
  - Smooth and supple texture with natural ability to sweat and easy to wash;

- Duvalcare The fabric is breathable, which prevents the occurrence of moisture and odors;
- It can be washed to a maximum temperature of 60 ° C in delicate wash cycle. It can also be dry-cleaned;
- You can choose between three densities depending on body weight or user preferences:
  - Normal: Latex with density 60kg/ m<sup>3</sup> for users weighing up to 70kg or preference for low consistency mattresses. Memory foam with a density of 53,0kg/ m<sup>3</sup>.
  - Firm: Latex with a density of 70kg / m<sup>3</sup>, for users with a body weight up to 85kg, or preference for medium consistency mattresses. Memory foam with a density of 53,0kg / m<sup>3</sup>.
  - Extra firm: Latex with a density of 80kg / m<sup>3</sup>, for users with body weight above 85kg or preference for high consistency mattresses. Memory foam with a density of 53,0kg / m<sup>3</sup>.
- Soft light dimensions: 190x90x15 / 18cm;
- Soft venus dimensions: 190x90x17cm;
- Notes:
  - References mention mattresses for single beds (up to 190 x 90 cm);
  - There is also for couples (up to 200 x 150 cm);
  - Mattresses made by measure are available (by budget).

Code	Description	Selling unit
61 CL CD SU15 S	MATTRESS SOFT LIGHT 190X90X15 CM	CX. 1
61 CL CD SU18 S	MATTRESS SOFT LIGHT 190X90X18 CM	CX. 1
61 CL CD VE N	MATTRESS SOFT VENUS 190X90X17 CM	CX. 1