

## AQUACEL AG FOAM



### MAIN SPECIFICATIONS

- Aquacel Foam dressings offer comfort, simplicity and the benefits of contact healing of a hydrofiber layer.
- Reduces pain while in contact with wound.
- It does not cause injury to the wound when changing the bandage.
- It can be used to treat people with skin lesions, including pressure ulcers.
- Greater absorption and retention of fluids.
- The top layer is impermeable. This allows the control of moisture of the wound, avoiding viral and bacterial contamination. Allows the patient to bathe or soak with the dressing.
- The foam layer is soft, leaving the dressing more comfortable for the patient, in addition to absorbing and retaining excess fluid in the wound.
- A layer of hydrophobic with ionic silver, which on contact with the wound provides the perfect environment for healing.
- Available with and without adhesive. The presentation with adhesive, this consists of soft silicone, with good adhesion and easy removal. Avoid skin irritation or sensitivity.
- Layer with Silver, aiding in the fight against bacteria and facilitating healing.

Code

Description

Selling unit

01101022	AQUACEL AG FOAM ADHERENT 10X10 CM	CX. 10
01101010	AQUACEL AG FOAM ADHERENT 12,5X 12,5 CM	CX. 10
01101011	AQUACEL AG FOAM NON ADHERENT 10X10 CM	CX. 10