

AQUACEL FOAM



MAIN SPECIFICATIONS

- Foam dressing that offers comfort, plus benefits for healing, compound with multiple layers.
- Protective upper layer that provides evaporation of excess moisture and acts as a barrier impervious to viruses and bacteria.
- Smooth foam layer that increases patient comfort and absorbs excess fluid.
- Hydrofiber layer that gels on contact with wound exudate to maintain a moist environment. Blocks exudate from the wound, including its harmful components. The vertical absorption avoids the lateral expansion of fluid, reducing the risk of maceration.
- It adapts to the bed of the wound, minimizing the spaces where the bacteria can develop.
- Presentations with and without adhesive or border.
- In the presentation with adhesive: Smooth adhering silicone collar designed to adhere to the perilesional skin, not adhering to the wound bed.
- Adhesive that respects the perilesional skin, allowing an easy application and removal. Low potential for skin irritation or contact dermatitis.

Code	Description	Selling unit
01101014	AQUACEL FOAM ADHERENT 10X 10 CM	CX. 10

FACTORMED
Produtos Médicos & Hospitalares