

CUTICELL CLASSIC PARAFFIN GAUZE



MAIN SPECIFICATIONS

- The paraffin-impregnated dressing.
- Cuticell Classic is the economical solution when treating exudating wounds. The low-adherent dressing is made of open weave cotton gauze impregnated with soft paraffin – to facilitate dressing change and to support a moist environment for the wound.
- Impregnated with soft paraffin
- Open weave 100% cotton fabric with interlocking threads
- Also available in handy dispenser box Classic offers reliable treatment quality.
- Low-adherent and soothing
- Facilitates dressing change
- Helps maintain a moist environment to support wound healing
- For increased stability when cut to size
- Does not fray
- Allows exudate to drain freely into an absorbent secondary dressing
- Prevents occlusive effects
- Economical treatment
-

Benefit of simple handling

- Reducing packaging waste
- Indications: For exudating wounds such as First and second-degree burns, Lacerations and abrasions, Leg ulcers, Donor and recipient skin graft sites

Code	Description	Selling unit
01254011	CUTICELL CLASSIC 5X5 CM	CX. 50
01254012	CUTICELL CLASSIC 10x10 CM	CX. 100