

LEUKOPLAST S LF PERFORATED ADHESIVE



MAIN SPECIFICATIONS

- Perforated roll adhesive suitable for fixing the skin material (eg. catheters, needles, tubes, probes, dressings, etc.).
- Composition synthetic rubber adhesive and adhesive synthetic resins, zinc oxide.
- Supported on paper.
- White colour.
- Good adhesion.
- Well tolerated by the skin allows the skin to breathe.

Code	Description	Selling unit
3254005	LEUKOPLAST S LF PERFORATED ADHESIVE WHITE 5 M X 10 CM	CX. 1
3254006	LEUKOPLAST S LF PERFORATED ADHESIVE WHITE 5 M X 18 CM	CX. 1