

# TENSOPLAST SPORT ELASTIC AND ADHESIVE BANDAGE



## MAIN SPECIFICATIONS

- Porous elastic adhesive bandage with protection band along the entire band
- Extensible lengthwise by 50%
- High porosity that allows perspiration, minimizes the risk of maceration
- Hypoallergenic
- Water resistant, water repellent
- Latex free

Code	Description	Selling unit
04254010	TENSOPLAST SPORT BANDAGE 2,5 M X 10 CM	CX. 1
04254008	TENSOPLAST SPORT BANDAGE 2,5 M X 6 CM	CX. 1
04254009	TENSOPLAST SPORT BANDAGE 2,5 M X 8 CM	CX. 1