

# PREVIZINC E BANDAGE WITH OXID ZINC PASTE



## MAIN SPECIFICATIONS

- Elastic bandage of zinc oxide paste with a longitudinal elasticity of approximately 80%, composed of 46% viscose and 54% polyamide with 15% zinc oxide paste.
- Texture: Double weave cloth soaked in zinc oxide paste.
- DTEX polyamide fold 78x29.5 yarn / cm Viscose weft NM 6 7 yarn / cm; 27 g / m<sup>2</sup> fabric when stretched; Zinc oxide paste 180 g / m<sup>2</sup> when stretched.
- Ligatures individually wrapped in aluminum foil bag according to product specification.
- Previzinc is a bandage impregnated with a zinc oxide composite paste, which is suitable for the treatment of thrombophlebitis, chronic venous insufficiency, bumps, fractures and traumas, dislocations of contusions and distortions.

Code	Description	Selling unit
4335000	PREVIZINC BANDAGE 10 M X10 CM	CX. 8
4335001	PREVIZINC E BANDAGE 5 M X10 CM	CX. 10
4335002	PREVIZINC E BANDAGE 7 M X10 CM	CX. 10

FACTORMED  
Produtos Médicos & Hospitalares