

HYDROSORB HYDROGEL DRESSING



MAIN SPECIFICATIONS

- Hydrocellular gel dressing with semipermeable top layer, but impermeable to germs; Transparent.
- It offers a balanced humidity environment, due to high concentration of water, about 60%.
- Promotes granulation and epithelialization.
- Excellent biocompatibility.
- Good cushioning effect.
- Can be removed without irritating the skin.
- Available with adhesive film on the lips.
- The gel dressing has the self-adhesive effect; However for permanent fixation, additional safety with hypoallergenic adhesive tapes or bandage is recommended.
- Hydrosorb comfort has edges with impermeable hypoallergenic adhesive film for secure attachment.
- For the wet treatment of wounds, especially for the treatment of non-infectious wounds clinically during the granulation and epithelization stage, in leg ulcers or decubits; For the treatment of burns up to grade 2; Promotes re-epitalization of the skin, to cover areas where skin excerpts have been removed, after the bleeding has stopped; Has a detachable sheet where the size of the wound can be written to be documented.

Code	Description	Selling unit
9008531 (H)	HYDROSORB 5x7,5 CM	CX. 5
9008541 (H)	HYDROSORB 10X10 CM	CX. 5
9008551 (H)	HYDROSORB 20X20 CM	CX. 3
9007021 (H)	HYDROSORB COMFORT 4,5x6,5 CM	CX. 5
9007061 (H)	HYDROSORB COMFORT 7,5x10 CM	CX. 5
1046059	HYDROSORB COMFORT 12,5X12,5 CM	CX. 5