

HYDROCOLL HYDROCOLLOID DRESSING



MAIN SPECIFICATIONS

- Hydrocoll is a self-adhesive absorbent hydrocolloid dressing with a semipermeable cover layer, but impermeable to bacteria and liquids.
- The hydrocolloid turns into a gel when it comes into contact with the wound exudate and forms a moist environment that stimulates the effect of granulation and epithelization.
- Available in different sizes and forms for every need.
- For the humid treatment of clinically uninfected wounds, wounds with mild to moderate exudate, especially in wounds prone to poor treatment and long duration for granulation formation, for example, leg ulcers or decoctions; For burns of degree 2; Hydrocoll sacral is used in case of decumbent ulcer in sacral region, Hydrocoll concave is used in elbows or heel.

Code	Description	Selling unit
1046001	HYDROCOLL 10X10 CM	CX. 10
1046053	HYDROCOLL 15X15 CM	CX. 5
1046002	HYDROCOLL 7,5X7,5 CM	CX. 10
1046056	HYDROCOLL CLASSIC 15X15 CM	CX. 10
1046032	HYDROCOLL FINO 10X10 CM	CX. 10

1046031	HYDROCOLL THIN 7,5X7,5 CM	CX. 10
1046060	HYDROCOLL STANDARD SACRAL 12X18 CM	CX. 5
9007492 (H)	HYDROCOLL 20X20 CM	CX. 5
9007402 (H)	HYDROCOLL 5X5 CM	CX. 10
9009381 (H)	HYDROCOLL CLASSIC 10X10 CM	CX. 10
9009401 (H)	HYDROCOLL CLASSIC 20X20 CM	CX. 10
9007562 (H)	HYDROCOLL STANDARD CONCAVE 8X12 CM	CX. 10
9007602 (H)	HYDROCOLL THIN 15X15 CM	CX. 5