

# HYDROTAC FOAM DRESSING



## MAIN SPECIFICATIONS

- A foam dressing which is coated with a net-like hydrogel layer on the side facing towards the wound and a gas permeable, waterproof and bacteria-resistant outer layer made of polyurethane; the absorbent foam material quickly absorbs wound exudate without drying out the wound; in relatively dry wounds the hydrogel releases moisture and promotes the formation and maintenance of a moist wound environment, whereby the granulation and epithelisation are promoted; the wound dressing is soft and conformable, provides a good cushioning effect and can be removed almost painlessly; also available with an adhesive border.
- HydroTac comfort has an adhesive border all around, which is coated with hypoallergenic polyacrylate adhesive and holds the dressing securely in place.
- For the treatment of lightly to medium exuding wounds, during the granulation and epithelization phases, e.g. for leg ulcers, diabetic foot syndrome, decubital ulcer; for the treatment of burns up to degree 2a and skin graft donor sites after adequate bleeding control.

Code	Description	Selling unit
6858498 (H)	HYDROTAC 6 CM Ø	CX. 10
6858329 (H)	HYDROTAC 10X10 CM	CX. 10
6858378 (H)	HYDROTAC 12,5X12,5 CM	CX. 10

6858409 (H)	HYDROTAC 15X15 CM	CX. 10
6858418 (H)	HYDROTAC 10X20 CM	CX. 10
6858438 (H)	HYDROTAC 20X20 CM	CX. 10
6858258 (H)	HYDROTAC CONCAVE 18X18,5 CM	CX. 10
6858108 (H)	HYDROTAC COMFORT 8X8 CM	CX. 10
6858158 (H)	HYDROTAC COMFORT 12,5X12,5 CM	CX. 10
6858188 (H)	HYDROTAC COMFORT 15X15 CM	CX. 10
6858138 (H)	HYDROTAC COMFORT 10X20 CM	CX. 10
6858218 (H)	HYDROTAC COMFORT 15X20 CM	CX. 10
6858128 (H)	HYDROTAC COMFORT 20X20 CM	CX. 10
6858278 (H)	HYDROTAC COMFORT SACRAL 18X18 CM	CX. 10
6858508 (H)	HYDROTAC COMFORT SACRAL 22 X 22 CM	CX. 10