

ASKINA ABSORB+ ABSORBENT DRESSING



MAIN SPECIFICATIONS

- Askina Absorb is a sterile dressing consisting of an inner layer of polymer fibers with strong water retention properties. The layer is placed in a non-woven polyester / polyethylene bag.
- In the presence of exudate, the polymer fibers form a moist gel which helps to maintain the wet wound environment by providing the natural healing conditions.
- Indications: Wounds with moderate to high amount of exudate.

Code	Description	Selling unit
3331010	ASKINA ABSORB+ 10X10 CM	CX. 10
3331210	ASKINA ABSORB + 10X20 CM	CX. 10
3332005	ASKINA ABSORB + 20X20 CM	CX. 10