

# ASKINA CALGITROL PASTE SILVER AND ALGINATE



## MAIN SPECIFICATIONS

- Askina Calgitrol Ag Paste is compounded with silver alginate matrix for medical use. High conformability allows for closer contact between an active matrix of silver and rubber and wound rust, in fact important when there are difficult to access wounds, such as tunnels and sinuses, common in severe burns and diabetic foot ulcers.
- Indications: slightly exudative wounds in granulation, primary alleviation in pressure ulcers (Stage II to IV).
- Healing enhancement - Accelerated healing: has adaptation to the contours of the lesion, facilitating the absorption of exudate and accelerating a cellular proliferation.
- Go to other products of the same sex - The cannula allows the product application in locations, moisturizing and healing in depth.
- Askina Calgitrol Highly compatible paste for direct contact with the wound bed; 100% silver ions, active against a broad spectrum of microorganisms, inc. MRSA; No need for activation; Continuous and controlled release for protection; Prolonged anti-bacterial action; Maintains hydration of the wound; Easy to remove by irrigation; Leaves no stains.

Code

Description

Selling unit

01001002

ASKINA CALGITROL PASTE 15 GR

CX. 10

---

6241505

ASKINA CALGITROL PASTE 15 GR

CX. 5

---