

ASKINA DERM STERILE FILM



MAIN SPECIFICATIONS

- Polyurethane film, coated with acrylic adhesive, elastic, sterile and smooth for skin, self-adhesive easy to apply. Capable of keeping the moist environment ideal for healing.
- Features: Innovative application system; High moisture vapor transmission rate (MVTR); Low adhesion adhesive, avoiding aggressive removal; Transparent, allows visual inspection of the wound .; Impervious to liquids and bacteria.
- Indication: Askina Derm is recommended for the management of non-exudative wounds, such as: Stage I ulcers; Pressure ulcers; Donor areas; Small abrasions; First degree burns; Surgical wounds; Lacerations; Protection against skin breakage due to friction or continuous exposure to moisture; Puncture site protection.
- Contraindication: Does not present.

| Code | Description | Selling unit |
|--------|----------------------|--------------|
| F72036 | ASKINA DERM 10X12CM | CX. 50 |
| F72037 | ASKINA DERM 10X20CM | CX. 10 |
| F72038 | ASKINA DERM 15X20 CM | CX. 10 |
| F72033 | ASKINA DERM 6X7CM | CX. 100 |