

# ASKINA DRESSIL SILICONE AND FOAM DRESSING



## MAIN SPECIFICATIONS

- Askina Dressil is a self adherent foam dressing with soft silicone adhesive on one side and a vapour permeable waterproof film on the other. It combines the absorption capacity of the foam with the soft adhesion of the silicone contact layer.
- Gentle and secure adherence.
- Minimized trauma during dressing changes.
- Significant absorption capacity.
- Easy to use.
- Good conformability on difficult-to-dress areas and skin contours.
- Repositionable, can be cut to shape.
- Improves patient comfort.
- Hypoallergenic.
- Indication: Pressure ulcers, venous leg ulcers, arterial leg ulcers, diabetic foot ulcers, traumatic wounds, surgical wounds, 1st and 2nd degree burns, donor sites and skin tears.
- Available with and without border.

Code	Description	Selling unit
5295710	ASKINA DRESSIL 5X7CM	CX. 10
5291010	ASKINA DRESSIL 10X10CM	CX. 10
5291210	ASKINA DRESSIL 10X20CM	CX. 10
5291510	ASKINA DRESSIL 15X15CM	CX. 10
5292005	ASKINA DRESSIL 20X20CM	CX. 5
5391010	ASKINA DRESSIL BORDER 10X10CM	CX. 10
5391503	ASKINA DRESSIL BORDER 15X15CM	CX. 10
5395210	ASKINA DRESSIL BORDER 15X20CM	CX. 10
5392005	ASKINA DRESSIL BORDER 20X20CM	CX. 5
5397510	ASKINA DRESSIL BORDER 7.5X7.5CM	CX. 10
5491605	ASKINA DRESSIL SACRUM 16X17.5CM	CX. 5