

## ASKINA HEEL FOAM DRESSING



### MAIN SPECIFICATIONS

- Askina Heel is an anatomically shaped two-layered, non-adherent foam dressing, which includes a foam strap and two self adhesive hooks that allow the dressing to remain in place even during movement.
- Indications: Management of moderate to heavily exuding, partial to full thickness wounds.
- Features and Benefits: Anatomically shaped; Easy to apply, providing protection to the malleolus; Non adhesive foam; Less trauma to surrounding skin, painless dressing changes; Self adhesive foam strap; Offers alternative fixation to hold the dressing in place.

Code	Description	Selling unit
7240110	ASKINA HEEL	CX. 10