

# ASKINA SOFT IV ADHESIVE



## MAIN SPECIFICATIONS

- Askina soft I.V. is a sterile and hypoallergenic dressing for cannula fixation, with a non-adherent and integrated absorbent pad
- Non-stick and absorbent pad:
- The puncture site remains dry
- Prevents grip at the puncture site
- Hypoallergenic acrylic adhesive
- Breathable fabric
- Cannula to tighten the cannula tube
- Extra cushion for extra comfort
- Indications: Fixation dressing for I.V.

Code	Description	Selling unit
9086919	ASKINA SOFT I.V. 8 X 6 CM	CX. 50