

# ASKINA SKIN FREEZE CRYOGENIC SPRAY



## MAIN SPECIFICATIONS

- Askina skin freeze is a cryogenic product for the treatment of warts and benign skin lesions.
- It has as an indication the removal of warts comfortably, in a safe and effective way.
- Warts treatment agent, consisting of a mixture of dimethyl ether and propane, which do not harm the ozone layer.
- It is the ideal cryotherapy treatment for benign skin lesions, combining the effectiveness of traditional cryosurgery with the advantages of an aerosol system application. Easy and efficient, without investing in equipment.
- Contains 50 disposable applicators of 5mm or 50 disposable applicators of 2mm, depending on the model, allowing a more accurate application and avoids the risk of cross contamination. Reaches a temperature of -55 ° C.
- Advantages: Easy to use; Used for various treatments; Minimum post-treatment complications; Compact container; There are no evaporative losses; Easy storage and transportation; Contains 2 applicators of different sizes.
- The effective treatment of: common warts, seborrheic fibroids, plantar warts, flat warts condyloma acuminata; Actinic keratoses.

Code

Description

Selling unit

9380809	ASKINA SKIN FREEZE 150 ML 5 MM	CX. 1
---------	--------------------------------	-------

---

9380701	ASKINA SKIN FREEZE SPRAY 150 ML SMALL 2 MM	CX. 1
---------	--	-------

---