

# AQUAFLATE GLYCERIN CATHETER SILICONE 2 WAYS



## MAIN SPECIFICATIONS

- Foley cushion AQUAFLATE Glycerine BRILLANT in silicone, transparent, with pre-filled syringe for balloon filler.
- Pre-filled syringe in sterile carton (10 ml of sterile 10% glycerol solution).
- Valve for Luer and Luer-lock syringe tips.
- Radiopaque tip and contrast line along the harness.
- Sterile, single use.
- Without latex.

| Code      | Description                                       | Selling unit |
|-----------|---------------------------------------------------|--------------|
| 178000-12 | CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 12 | CX. 5        |
| 178000-14 | CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 14 | CX. 5        |
| 178000-16 | CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 16 | CX. 5        |
| 178000-18 | CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 18 | CX. 5        |
| 178000-20 | CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 20 | CX. 5        |

|           |                                                   |       |
|-----------|---------------------------------------------------|-------|
| 178000-22 | CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 22 | CX. 5 |
|-----------|---------------------------------------------------|-------|

---

|           |                                                   |       |
|-----------|---------------------------------------------------|-------|
| 178000-24 | CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 24 | CX. 5 |
|-----------|---------------------------------------------------|-------|

---