

AQUAFLATE GLYCERIN CATHETER SILICONE 2 WAYS



MAIN SPECIFICATIONS

- Foley cushion AQUAFLATE Glycerine BRILLANT in silicone, transparent, with pre-filled syringe for balloon filler.
- Pre-filled syringe in sterile carton (10 ml of sterile 10% glycerol solution).
- Valve for Luer and Luer-lock syringe tips.
- Radiopaque tip and contrast line along the harness.
- Sterile, single use.
- Without latex.

Code	Description	Selling unit
178000-12	CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 12	CX. 5
178000-14	CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 14	CX. 5
178000-16	CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 16	CX. 5
178000-18	CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 18	CX. 5
178000-20	CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 20	CX. 5

178000-22	CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 22	CX. 5
-----------	---	-------

178000-24	CATHETER AQUAFLATE GLYCERIN SILICONE 2 WAYS CH 24	CX. 5
-----------	---	-------
