

URINARY CATHETER NELATION STERILE FLOCATH HYDROGEL



MAIN SPECIFICATIONS

- The basic Flocath Hydrogel offers all the advantages of the Flocath hydrophilic technology for a high standard safe and comfortable catheterisation.
- Approved safety.
- The Flocath catheter made of DEHP-free PVC with vertically cut and softly rounded eyes preserves the sensitive urethral mucosa.
- The silky hydrophilic coating of the catheter, based on PVP (polyvinyl pyrrolidone), provides a fast and homogeneous hydration, a high biocompatibility, and is resistant to abrasion. Its superior lubricant qualities allow it to slide easily into position.
- Comfortable application.
- The special Flocath coating offers significantly lower friction than conventional coatings.
- extended dry-out time up to 10-15 minutes allows a relaxed catheterisation.
- Clean – no gel residue left in the urethra, bladder or on the hands.
- Easy to handle – the packaging hole enables the pack to be opened easily, even with limited manual dexterity.

- Convenient – the self-adhesive patch and the packaging hole allow the pack to be fixed in place or hung up where convenient while preparing the catheterisation procedure with both hands free.
- Additional – to activate the hydrophilic coating the catheter has to be moistened by, preferably, sterile water or saline solution.
- Individual – Flocath Hydrogel is available in different lengths, diameters (sizes) and catheter tips.
- Expandable – if a drainage bag is needed, it can be connected with the funnel at the end of the catheter.
- Gentle on the urethra and the urethral mucosa thanks to the Flocath catheter.
- Smooth and comfortable application due to an evenly applied silky and biocompatible hydrophilic coating.
- Stressless catheterisation through an extended dry-out time of the coating.
- Safe, clean and economic.
- Length 20 cm for women and 40 cm for men. Available in several sizes.

Code	Description	Selling unit
85 16 41 - 8	FLOCATH HYDROGEL 40 CM MASCULINE CH 8	CX. 30
85 16 41 - 10	FLOCATH HYDROGEL 40 CM MASCULINE CH 10	CX. 30
85 16 41 - 12	FLOCATH HYDROGEL 40 CM MASCULINE CH 12	CX. 30
85 16 41 - 14	FLOCATH HYDROGEL 40 CM MASCULINE CH 14	CX. 30
85 16 41 - 16	FLOCATH HYDROGEL 40 CM MASCULINE CH 16	CX. 30
85 16 41 - 18	FLOCATH HYDROGEL 40 CM MASCULINE CH 18	CX. 30
85 16 21 - 8	FLOCATH HYDROGEL 20 CM FEMININE CH 8	CX. 30
85 16 21 - 10	FLOCATH HYDROGEL 20 CM FEMININE CH 10	CX. 30
85 16 21 - 12	FLOCATH HYDROGEL 20 CM FEMININE CH 12	CX. 30
85 16 21 - 14	FLOCATH HYDROGEL 20 CM FEMININE CH 14	CX. 30
85 16 21 - 16	FLOCATH HYDROGEL 20 CM FEMININE CH 16	CX. 30