

CATHETER URINARU NELATON STERILE

FLOCATH QUICK



MAIN SPECIFICATIONS

- The FloCath Quick Hydrophilic Intermittent Catheter is an "all-in-one" intermittent catheter providing the ultimate convenience.
- The integrated package contains a FloCath Hydrophilic Catheter with our patented handling sheath and a packet of sterile 0.9% saline.
- Before opening, break the sterile saline pouch to hydrate the catheter.
- The lubricious coating allows for a comfortable insertion.
- The handling sheath provides an easy grip, allowing for touchless insertion.
- The sheath can also be pushed back and used as an extension to the catheter.
- Length of 20 cm for women and 40 cm for men. Available in several sizes.

Code	Description	Selling unit
85 12 41 - 8	FLOCATH QUICK 40 CM MASCULINO CH 8	CX. 30
85 12 41 - 10	FLOCATH QUICK 40 CM MASCULINE CH 10	CX. 30
85 12 41 - 12	FLOCATH QUICK 40 CM MASCULINE CH 12	CX. 30

85 12 41 - 14	FLOCATH QUICK 40 CM MASCULINE CH 14	CX. 30
85 12 41 - 16	FLOCATH QUICK 40 CM MASCULINE CH 16	CX. 30
85 12 41 - 18	FLOCATH QUICK 40 CM MASCULINE CH 18	CX. 30
85 12 41 - 20	FLOCATH QUICK 40 CM MASCULINE CH 20	CX. 30
85 12 21 - 8	FLOCATH QUICK 20 CM FEMININE CH 8	CX. 30
85 12 21 - 10	FLOCATH QUICK 20 CM FEMININE CH 10	CX. 30
85 12 21 - 12	FLOCATH QUICK 20 CM FEMININE CH 12	CX. 30
85 12 21 - 14	FLOCATH QUICK 20 CM FEMININE CH 14	CX. 30
85 12 21 - 16	FLOCATH QUICK 20 CM FEMININO CH 16	CX. 30
85 12 21 - 18	FLOCATH QUICK 20 CM FEMININO CH 18	CX. 30
85 12 21 - 20	FLOCATH QUICK 20 CM FEMININO CH 20	CX. 30