

# SURGICAL GLOVES LATEX WITHOUT POWDER GYNOGLOVE



## MAIN SPECIFICATIONS

- Surgical gloves with a lengthened cuff.
- Glove type: surgical, sterile, with a lengthened cuff
- Raw material: natural rubber latex
- Powdering agent: none
- Quantity in a packaging: 1 pair, 25 pairs, 100 pairs
- Surface: texture on the prehensile side of the hand
- Cuff finishing: rolled edge
- Shape: anatomical, different for the right and left hand
- Sizes: 6.5 (S), 7.5 (M), 8.5 (L)

Code	Description	Selling unit
-	SURGICAL GLOVES LATEX WITHOUT POWDER GYNOLOVE NR 6,5	CX. 50
--	SURGICAL GLOVES LATEX WITHOUT POWDER GYNOGLOVE NR 7.5	CX. 50

---

SURGICAL GLOVES LATEX WITHOUT POWDER GYNOGLOVE  
NR 8,5

CX. 50

---

FACTORMED  
Produtos Médicos & Hospitalares