

# SURGICAL GLOVES LATEX WITH POWDER

## INTOUCH



### MAIN SPECIFICATIONS

- Surgical pre-powdered latex gloves with a straight cuff
- glove type: surgical, sterile
- raw material: natural rubber latex
- powdering agent: corn starch (meal)
- quantity in a packaging: 1 pair, 50 pairs, 200 pairs
- surface: micro-textured
- cuff finishing: straight, without rolled finishing
- shape: anatomical, different for the right and left hand
- sizes: 6.0 - 9.0

Code	Description	Selling unit
-	SURGICAL GLOVES LATEX WITH POWDER INTOUCH NR 6	CX. 50
--	SURGICAL GLOVES LATEX WITH POWDER INTOUCH NR 6,5	CX. 50
---	SURGICAL GLOVES LATEX WITH POWDER INTOUCH NR 7	CX. 50

----	SURGICAL GLOVES LATEX WITH POWDER INTOUCH NR 7,5	CX. 50
-----	SURGICAL GLOVES LATEX WITH POWDER INTOUCH NR 8	CX. 50
-----	SURGICAL GLOVES LATEX WITH POWDER INTOUCH NR 8,5	CX. 50
-----	SURGICAL GLOVES LATEX WITH POWDER INTOUCH NR 9	CX. 50