

SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT



MAIN SPECIFICATIONS

- Dark-blue gloves are excellent for the double gloving system and high-risk procedures.
- Glove type: surgical, sterile
- Raw material: natural rubber latex
- Powdering agent: none - with an internal latex-free layer
- Quantity in a packaging: 1 pair, 50 pairs, 200 pairs
- Surface: micro-textured
- Cuff finishing: rolled edge
- Shape: anatomical, different for the right and left hand
- Sizes: 6.0 - 9.0

Code	Description	Selling unit
-	SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT NR 6	CX. 50
--	SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT NR 6,5	CX. 50

---	SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT NR 7	CX. 50
----	SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT NR 7,5	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT NR 8	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT NR 8,5	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER INTOUCH SPOT NR 9	CX. 50