

SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx



MAIN SPECIFICATIONS

- Brown micro-surgical gloves.
- glove type: surgical, sterile, with very thin walls
- raw material: natural rubber latex
- powdering agent: none
- quantity in a packaging: 1 pair, 50 pairs, 300 pairs
- surface: textured
- cuff finishing: rolled edge
- shape: anatomical, different for the right and left hand
- sizes: 6.0 - 9.0

Code	Description	Selling unit
-	SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx NR 6	CX. 50
--	SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx NR 6,5	CX. 50

---	SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx NR 7	CX. 50
----	SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx NR 7,5	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx NR 8	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx NR 8,5	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER MICROTEx NR 9	CX. 50