

SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG



MAIN SPECIFICATIONS

- Brown orthopaedic gloves.
- Glove type: surgical, sterile, orthopaedic, with thicker walls
- Raw material: natural rubber latex
- Powdering agent: none
- Quantity in a packaging: 1 pair, 50 pairs, 300 pairs
- Surface: textured
- Cuff finishing: rolled edge
- Shape: anatomical, different for the right and left hand
- Sizes: 6.0 - 9.0

Code	Description	Selling unit
-	SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG NR 6	CX. 50
--	SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG NR 6,5	CX. 50

---	SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG NR 7	CX. 50
----	SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG NR 7,5	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG NR 8	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG NR 8,5	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER ORTHOPEG NR 9	CX. 50