

# SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE



## MAIN SPECIFICATIONS

- Type of the glove: Powder free, sterile surgical glove
- Material Natural rubber latex
- Donning powder: none
- Colour: Creamy
- Shape: Anatomic, hand specific
- Cuff: Beaded
- External surface: Textured
- Sizes: 6.0 - 9.0

Code	Description	Selling unit
-	SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE NR 6	CX. 50
--	SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE NR 6,5	CX. 50
---	SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE NR 7	CX. 50
----	SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE NR	CX. 50

7.5

-----	SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE NR 8	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE NR 8,5	CX. 50
-----	SURGICAL GLOVES LATEX WITHOUT POWDER HANDSAFE NR 9	CX. 50