

# SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC NEOPRENE



## MAIN SPECIFICATIONS

- Surgical neoprene gloves
- Glove type: surgical, sterile
- Raw material: neoprene
- Powdering agent: none
- Quantity in a packaging: 1 pair, 50 pairs, 400 pairs
- Surface: textured
- Cuff finishing: rolled edge
- Shape: anatomical, different for the right and left hand
- Sizes: 6.0 - 9.0

| Code | Description  | Selling unit |
|------|--|--------------|
| -    | SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC NEOPRENE NR 6   | CX. 50       |
| --   | SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC NEOPRENE NR 6,5 | CX. 50       |

|       |   |        |
|-------|---|--------|
| ---   | SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC<br>NEOPRENE NR 7   | CX. 50 |
| ----  | SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC<br>NEOPRENE NR 7,5 | CX. 50 |
| ----- | SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC<br>NEOPRENE NR 8   | CX. 50 |
| ----- | SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC<br>NEOPRENE NR 8,5 | CX. 50 |
| ----- | SURGICAL GLOVES NEOPRENE WITHOUT POWDER SYNTEC<br>NEOPRENE NR 9   | CX. 50 |