

EXAM GLOVES LATEX WITH POWDER

COMFORT



MAIN SPECIFICATIONS

- Pre-powdered latex gloves, diagnostic and protective
- glove type: diagnostic, protective, non-sterile
- Gaw material: natural rubber latex
- Powdering agent: corn starch (meal)
- Quantity in an individual packaging: 100 pieces
- Surface: smooth
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	EXAM GLOVES LATEX WITH POWDER COMFORT XS	CX. 100
--	EXAM GLOVES LATEX WITH POWDER COMFORT S	CX. 100
---	EXAM GLOVES LATEX WITH POWDER COMFORT M	CX. 100

----	EXAM GLOVES LATEX WITH POWDER COMFORT L	CX. 100
----	EXAM GLOVES LATEX WITH POWDER COMFORT XL	CX. 100