

EXAM GLOVES NITRILE WITHOUT POWDER

NITRYLEX



MAIN SPECIFICATIONS

- Powder-free nitrile gloves recommended for persons allergic to latex
- Glove type: diagnostic, protective, non-sterile
- Raw material: nitrile
- Powdering agent: none
- Colour: blue/purple
- Quantity in an individual packaging: 100 pieces
- Surface: textured at finger tips
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX XS	CX. 100
--	EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX S	CX. 100

---	EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX M	CX. 100
----	EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX L	CX. 100
-----	EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX XL	CX. 100