

EXAM GLOVES NITRILE WITHOUT POWDER

NITRYLEX CHEMO



MAIN SPECIFICATIONS

- Powder-free nitrile gloves recommended for persons allergic to latex
- Glove type: diagnostic, protective, non-sterile
- Raw material: nitrile
- Powdering agent: none
- Colour: blue
- Quantity in an individual packaging: 100 pieces
- Surface: textured at finger tips
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

| Code | Description | Selling unit |
|------|--|--------------|
| - | EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX CHEMO XS | CX. 100 |

| | | |
|-------|---|---------|
| -- | EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX CHEMO S | CX. 100 |
| --- | EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX CHEMO M | CX. 100 |
| ---- | EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX CHEMO L | CX. 100 |
| ----- | EXAM GLOVES NITRILE WITHOUT POWDER NITRYLEX CHEMO XL | CX. 100 |