

# EXAM GLOVES LATEX WITH POWDER

## TEXIDENT



### MAIN SPECIFICATIONS

- Pre-powdered latex gloves, diagnostic and protective, recommended for dentists
- glove type: diagnostic, protective, non-sterile
- raw material: natural rubber latex
- powdering agent: corn starch (meal)
- quantity in an individual packaging: 100 pieces
- surface: micro-textured
- cuff finishing: evenly rolled edge
- shape: universal, fitting the right and left hand
- sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	EXAM GLOVES LATEX WITH POWDER TEXIDENT XS	CX. 100
--	EXAM GLOVES LATEX WITH POWDER TEXIDENT S	CX. 100
---	EXAM GLOVES LATEX WITHO POWDER TEXIDENT M	CX. 100

----	EXAM GLOVES LATEX WITH POWDER TEXIDENT L	CX. 100
----	EXAM GLOVES LATEX WITH POWDER TEXIDENT XL	CX. 100