

# EXAM GLOVES VINYL WITH POWDER

## VINYLEX BLUE



### MAIN SPECIFICATIONS

- Blue vinyl gloves recommended for persons allergic to latex
- Glove type: diagnostic, protective, non-sterile
- Raw material: vinyl polychloride
- Powdering agent: corn starch (meal)
- Quantity in an individual packaging: 100 pieces
- Surface: smooth
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	EXAM GLOVES VINYL WITH POWDER VINYLEX BLUE XS	CX. 100
--	EXAM GLOVES VINYL WITH POWDER BLUE S	CX. 100
---	EXAM GLOVES VINYL WITH POWDER VINYLEX BLUE M	CX. 100

----	EXAM GLOVES VINYL WITH POWDER VINYLEX BLUE L	CX. 100
-----	EXAM GLOVES VINYL WITH POWDER VINYLEX BLUE XL	CX. 100