

DRUGWIPE 5 DRUG DETECTION TEST / SALIVA METHOD



MAIN SPECIFICATIONS

- Storage & Disposal

Always keep DrugWipe® 5 in environments between + 5 ° C and +25 ° C. Never expose the device to temperatures below 0 ° C (do not let it freeze) or above + 40 ° C. After use or in case of destruction, the devices must be recycled. Check the recycling procedures of your area of residence.

- Shipping & Handling

The transport time should not exceed 3 days. Keep the DrugWipe® 5 test at temperatures between + 5 ° C and + 35 ° C. When used on rainy or snowy days, take care that the device does not get any moisture. Testing should be done at temperatures between + 5 ° C and + 40 ° C. If the test is used at temperatures below + 5 ° C, it must be preheated (eg by hand) before opening.

- Quality control

DrugWipe®5 + has an integrated control system that checks the functionality of the test and its correct use.

Inquire in your company about the additional norms of use and the due control of quality.

- Cross Reactions

Substances in the DrugWipe® 5+ test that may cause cross-reactivity were analyzed in our laboratories. If requested, we may send the complete list of tested substances to our customers.

- The different products
- DrugWipe® Single Tests
- No. of article Type of test
- 001-CAN DrugWipe® Cannabis Sativa (hashish, marijuana)
- 002-COC DrugWipe® Cocaine
- 003-OPI DrugWipe® Opiates
- 004-AMP DrugWipe® Amphetamines (Ecstasy / Methamphetamine)
- 005-BEN DrugWipe® Benzodiazepines
- 006-MTD DrugWipe® Methadone
- DrugWipe® TWIN
- No. of article Type of test

0014-TWIN-C / A DrugWipe® TWIN Cannabis / Amphetamines

0023-TWIN-C / O DrugWipe® TWIN Cocaine / Opiates

0014N-TWIN-C / A DrugWipe® TWIN Cannabis / Amphetamine saliva test

0023N-TWIN-C / O DrugWipe® TWIN Cocaine / Opiate Saliva Test

- DrugWipe® 5

No. of article Type of test

008-DW-5 DrugWipe® 5

Cannabis, Amphetamines (Ecstasy / Methamphetamine), Cocaine, Opiates

009N-DW-5 + DrugWipe® 5+ saliva test

Cannabis, Amphetamines (Ecstasy / Methamphetamines), Cocaine, Opiates

Code	Description	Selling unit
20115000	TESTE P/DETECÇÃO DE DROGAS DRUGWIPE® 5+ MÉTODO SALIVA	CX.25