

HOUSEHOLD GLOVES LATEX WITH POWDER IDEALL



MAIN SPECIFICATIONS

- Pre-powdered latex gloves, diagnostic and protective
- Glove type: diagnostic, protective, non-sterile
- Raw material: natural rubber latex
- Powdering agent: corn starch (meal)
- Quantity in an individual packaging: 100 pieces
- Surface: smooth
- Cuff finishing: evenly rolled edge
- Shape: universal, fitting the right and left hand
- Sizes: XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)

Code	Description	Selling unit
-	HOUSEHOLD GLOVES LATEX WITH POWDER IDEALL XS	CX. 100
--	HOUSEHOLD GLOVES LATEX WITH POWDER IDEALL S	CX. 100
---	HOUSEHOLD GLOVES LATEX WITH POWDER IDEALL M	CX. 100

----	HOUSEHOLD GLOVES LATEX WITH POWDER IDEALL L	CX. 100
-----	HOUSEHOLD GLOVES LATEX WITH POWDER IDEALL XL	CX. 100