

SYRINGE INJ/LIGHT (WITHOUT NEEDLE, CATHETER CONE)



MAIN SPECIFICATIONS

- PP barrell with graduated scale
- PP plunger
- Latex free gasket
- Catheter cone
- Syringe type: standard syringes
- Safety mechanism: no
- Needle: no

Code	Description	Selling unit
-	SYRINGE INJ/LIGHT (WITHOUT NEEDLE, CATHETER CONE) 60 ML	CX. 20
--	SYRINGE INJ/LIGHT (WITHOUT NEEDLE, CATHETER CONE) 100 ML	CX. 25

FACTORMED
Produtos Médicos & Hospitalares